

# INTENSIVE CARE VEHICLE

Mobile, negative pressure, ICU  
unit & Professional Air  
Sterilization Systems



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗  
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗  
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
Global unified hotline: +86-532-88935593

**Haier Biomedical**  
Make Life Better

### Features

Mobile negative pressure ICU is a professional mobile ward equipped with the best nursing, treatment and rehabilitation services, which can provide an isolated treatment place for critically ill patients. The ward is equipped with air supply purification device and air exhaust purification device, and the virus killing efficiency is 99.99%.

The vehicle environment meets the requirements of hospital clean room cleanliness and air sanitation disinfection. The purification efficiency of 0.3um dust particles is above 99.97%, and the sterilization efficiency is above 99.9%. It has the functions of sterilization, anti-bacteria, constant temperature, constant humidity, constant pressure, heat recovery, high efficiency and energy saving, which can effectively guarantee the air quality and safety of the ward.



### Parameters

Overall dimensions of the vehicle (mm)	9990*2500*3840
Vehicle weight (kg)	18000
Chassis model	Shandeka ZZ5187XXYN711GF1K
Emission Standards	Chinese Standard VI、GB3847-2005
Carriage	steel structure + thermal insulation board
Engine displacement (ml)	6870



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
 Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
 No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
 URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
 Global unified hotline: +86-532-88935593

Haier Biomedical  
 Make Life Better