

MOBILE MEDICAL EXAMINATION UNIT

Vehicle heating, ventilation
environment, interior and
exterior lighting



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better

Features

Haier Biomedical according to the needs of X-ray machine, color ultrasound, ECG, biochemical, gynecological examination bed, traditional Chinese medicine diagnosis and treatment, eye, ear, nose and throat instruments and other equipment, to fully meet the needs of various physical examination, but also to show the overall image and comprehensive strength of the hospital.

Chassis optional: Shandeka, Dongfeng, Qingling, Benz. Interior style, layout according to customer needs tailored.



In vitro diagnostics

Auscultation

Color Doppler

Parameters

Overall dimensions of the vehicle (mm)	9940*2470*3780
Vehicle weight (kg)	13500
Wheelbase (mm)	5600
Chassis engine maximum power (kW)	155
Extended box width (mm)	650
generator set (kW)	6
Maximum speed (km/h)	89



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better