



DW-40L278

Model	DW-40L278
Climate Class	N
Temperature Range(°C)	-20~ -40
Power Supply(V/Hz)	220-240/50/60
Power(W)	310
Capacity(L/Cu.Ft)	278/9.8
Net/Gross Weight (approx)	115/135(kg) 253.5/297.6(lbs)
Interior Dimension (W*D*H)	520*435*1230(mm) 20.5*17.1*48.4(in)
Exterior Dimension (W*D*H)	736*660*1810(mm) 29.0*26.0*71.3(in)
Drawer Dimensions (W*D*H)	457*354*167(mm) 18.0*13.9*6.6(in)
Foot	Y
Drawers	6

Haier New Generation Upright Deep Freezer

Key Features

- Microprocessor control, adjustable temperature range : -20°C to -40°C
- LED digital display and increment at 0.1°C
- Permanently lubricated cooling fan for safety and longevity
- Rapid cooling with shelf evaporator
- Removable double outer door seal design, good sealing effect and energy saving

Ergonomic Design

- USB Data logging and temperature recorder are optional
- Padlock design and stainless steel cylinder
- Tracks to label drawers for product identification

Safety

- Multiple malfunction alarms include high/low temperature, sensor error, power failure, high ambient temperature, low battery
- Alarm Types: buzzer, flashing light, remote alarm



Reliable Performance



USB














Malfunction Alarms



Double Outer Door Seal



Product Comparison

Model	DW-40L278	DW-40L262
Temperature Range	-20~-40°C	-10~-40°C
Temperature Control	Microprocessor control and increment at 0.1°C 	Microprocessor control and increment at 1°C 
Evaporator	Evaporator made of embossed aluminum sheet 	Evaporator made of plastic dripped pipe 
Seal	Double Seal Structure 	Embedded Removable Seal Structure 
Lock	Door Lock+Padlock 	Padlock only 
Drawer	labled card slot 	NA 
Loading Capacity	50*6=300 bags(200ml/bag)	24*6+12=156 bags(200ml/bag)
Multiple Malfunction Alarm	High and low temperature, Sensor error, Door open, Power failure, Low battery, High ambient temperature	High and low temperature, Sensor error,
Alarm Type	Buzzer , flashing light and remote alarm	Buzzer and flashing light
Optional Parts	USB Interface Temperature Recorder 	NA

Notice

The new product will be launched in 30th Nov, 2014. We hope respectful partners to investigate and market, if any request about the product and targeted price, please directly contact Grace. Suggestions on the newsletter are also appreciated.

The upgraded will be more professional, more stable than the current ones. And if there are some questions about this, please contact with Grace Xiao (xiaozh@haier.com) directly. If you have suggestions on newsletter and have some ideas to share with global dealers, please also contact me. TKS.