

5G MOBILE STROKE UNIT

An ambulance equipped with a small CT and related equipment



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better

Features

5G Mobile Stroke Unit (5GMSU) combines Mobile Computed Tomography (CT), MCT and related inspection instruments are assembled in the ambulance. Under the 5G network, a treatment unit integrating neurological symptom examination, CT diagnosis and intravenous thrombolytic therapy is formed.

5G mobile stroke unit adopts mobile medical innovation technology, mainly including: small mobile CT, scanner, inspection and monitoring medical equipment with information transmission function, 5G network, consultation center, remote stroke treatment center and other mobile medical technology and auxiliary equipment.

Parameters

Overall dimensions of the vehicle (mm)	5990*2300*2900
Vehicle weight (kg)	6050
Chassis model	NJ1066EFC
engine model	F1CE8481K
Emission Standards	Chinese Standard VI
Maximum speed (km/h)	120
Seat (person)	7
CT	16 rows of mobile CT
medical cabin	Lead protection – protection level meets the requirements of GBZ264-2015 standard



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better